

Modellvarianten/ *Model versions:*

Beispielhafte Zusammensetzung des Artikelcodes/
Example for the formation of the article code:

RA631 K 3 2 Z 65 22 2 8
 ① ② ③ ④ ⑤ ⑥ ⑦ ⑧

① Ausführung/ Type	
	Code
Kostüm/ Costume	K
Weste/ Vest	W
Rock/ Skirt	R
Weste Victoria/ Vest Victoria	V

② Bleigleichwert-Front/ Lead Equivalent-Front	
	Code
0.25 mm Pb	2
0.35 mm Pb	3
0.50 mm Pb	5
0.25 mm Pb / 0.42 mm Pb (Levio)	6
0.25 mm Pb / 0.50 mm Pb (Levio)	7

③ Bleigleichwert-Rücken/ Lead Equivalent-Back	
	Code
0.25 mm Pb	2

④ Größen/ Sizes			
	Code		Code
XS	Z	XL	X
S	S	XXL	D
M	M	3XL	T
L	L		

⑤ Rocklängen/ Skirt Length	
	Code
65 cm	65
70 cm	70
75 cm	75
80 cm	80

⑤ Farbauswahl/ Colours			
	Code		Code
Curacao	22	Titan	29
Lagoon	25	Mandarin	30
Orchidee	27	Blue Admiral	31
Regatta	28	Baccara	32

⑥ Außenmaterial – Hülle/ Outer Material - Cover	
	Code
ComforTex HPMF	1
ComforTex HPMF Hybrid	2

⑦ Strahlenschutzmaterial – Innen/ X-Ray Protective Material - Inlay	
	Code
Novalite Plus; 0.25 mm Pb 0.35 mm Pb 0.42 mm Pb and 0.50 mm Pb	1
Leadfree 2.0; 0.25 mm Pb 0.35 mm Pb 0.42 mm Pb and 0.50 mm Pb	3

MAVIG GmbH
 X-Ray Protection and
 Medical Suspension Systems

Trade-Register Munich
 HRB-No. 76 403

Managing Directors:
 Christian Stoian
 Alexia Stoian-Lepère
 Philip Milbrdt

VAT-ID No.:
 DE 811268845