



EN

RP668

MAVIG



Models:	RP66852S	Children's X-Ray Protective Coat (SMALL)
	RP66852M	Children's X-Ray Protective Coat (MEDIUM)
	RP66852L	Children's X-Ray Protective Coat (LARGE)

Made in Germany

CE 0302

RP668 – Children's All-round Protective Coat

The cut combines ideal radiation protection with variable fit and high wearing comfort. The flexible material and generous width adjustment allow for a comfortable fit.

Children's X-ray protection, in particular, must be easy to put on and take off. So, you can adjust the RP668 children's coat quickly and securely using the hook-and-loop fasteners at the front. Children's X-ray protection is manufactured by MAVIG with the innovative Leadfree 2.0 radiation protection material to achieve the best possible weight ratio and thus impose only the most necessary physical strain on the children.

- **Protective Material** Leadfree 2.0
- **Lead Equivalent** Front protection 0.50 mm Pb / Back panels 0.25 mm Pb
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Plus (design Happy only), available in the colours:
Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or design: Happy
- **Size**

SMALL	Length: 60 cm, Hip girth: 68 - 75 cm for children 3 to 5 years old
MEDIUM	Length: 80 cm, Hip girth: 76 - 83 cm for children 6 to 8 years old
LARGE	Length: 84 cm, Hip girth: 84 - 95 cm for children 9 to 12 years old

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com